


# CURETTES


Curettes are precision instruments used for scraping or debriding bone and soft tissue in orthopedic procedures. Designed with sharp, scoop-shaped tips, they allow for controlled removal of unwanted tissue while preserving surrounding structures. Manufactured from surgical-grade stainless steel, our curettes offer excellent balance, durability, and tactile feedback for precise surgical work.



# BRUN CURETTES

## BRUN (SPRATT) CURETTES

Oval cup, hollow handle




CODE #	TIP SHAPE	SIZE	LENGTH
SC46511	Straight	#5/0	7" (18 cm)
SC46512	Straight	#4/0	7" (18 cm)
SC46513	Straight	#3/0	7" (18 cm)
SC46515	Straight	#2/0	7" (18 cm)
SC46517	Straight	#0	7" (18 cm)
SC46519	Straight	#1	7" (18 cm)
SC46521	Straight	#2	7" (18 cm)
SC46523	Straight	#3	7" (18 cm)
SC46525	Straight	#4	7" (18 cm)
SC46529	Straight	#5	7" (18 cm)
SC46531	Straight	#6	7" (18 cm)
SC46563	Straight	#3/0	9" (23 cm)
SC46564	Straight	#2/0	9" (23 cm)
SC46565	Straight	#0	9" (23 cm)
SC46566	Straight	#1	9" (23 cm)
SC46567	Straight	#2	9" (23 cm)
SC46568	Straight	#3	9" (23 cm)
SC46569	Straight	#4	9" (23 cm)
SC46570	Straight	#5	9" (23 cm)
SC46571	Straight	#6	9" (23 cm)

CODE #	TIP SHAPE	SIZE	LENGTH
SC49110	Angled 45 degrees	#3/0	9" (23 cm)
SC49111	Angled 45 degrees	#2/0	9" (23 cm)
SC49112	Angled 45 degrees	#0	9" (23 cm)
SC49113	Angled 45 degrees	#1	9" (23 cm)
SC49114	Angled 45 degrees	#2	9" (23 cm)
SC49115	Angled 45 degrees	#3	9" (23 cm)
SC49116	Angled 45 degrees	#4	9" (23 cm)
SC49117	Angled 45 degrees	#5	9" (23 cm)
SC49118	Angled 45 degrees	#6	9" (23 cm)

# BRUN (SPRATT) CURETTES

Round up, hollow handle



CODE #	TIP SHAPE	SIZE	LENGTH
SC49120	Straight	#3/0	6-3/4" (17 cm)
SC49121	Straight	#2/0	6-3/4" (17 cm)
SC49122	Straight	#0	6-3/4" (17 cm)
SC49123	Straight	#1	6-3/4" (17 cm)
SC49124	Straight	#2	6-3/4" (17 cm)
SC49125	Straight	#3	6-3/4" (17 cm)
SC49126	Straight	#4	6-3/4" (17 cm)

# BRUN CURETTES


8-7/8" (22.5 cm).



CODE #	DESCRIPTION
343500	Size 000
343501	Size 00
343501	Size 0
343503	Size 1
343504	Size 2
343505	Size 3
343506	Size 4

# BRUN CURETTES

Solid hexagonal handle, 9" (23 cm).



CODE #	TIP SHAPE	SIZE
SC46551	Straight	#3/0
SC46552	Straight	#2/0
SC46553	Straight	#0
SC46554	Straight	#1
SC46555	Straight	#2
SC46556	Straight	#3
SC46557	Straight	#4
SC46558	Straight	#5
SC46559	Straight	#6

CODE #	TIP SHAPE	SIZE
SC49150	Angled 45 degrees	#3/0
SC49151	Angled 45 degrees	#2/0
SC49152	Angled 45 degrees	#0
SC49153	Angled 45 degrees	#1
SC49154	Angled 45 degrees	#2
SC49155	Angled 45 degrees	#3
SC49156	Angled 45 degrees	#4
SC49157	Angled 45 degrees	#5
SC49158	Angled 45 degrees	#6

# COBB CURETTES

## COBB CURETTES, ALUMINUM HANDLE

All stainless steel working end and anodized aluminum, knurled handle, 11" (28 cm).



CODE #	TIP SHAPE	SIZE
SC46433	Angled	#3/0
SC46434	Angled	#2/0
SC46435	Angled	#0
SC46436	Angled	#1
SC46437	Angled	#2
SC46438	Angled	#3
SC46439	Angled	#4
SC46440	Angled	#5
SC46441	Angled	#6

CATALOG #	TIP SHAPE	SIZE
SC46422	Straight	#3/0
SC46423	Straight	#2/0
SC46424	Straight	#0
SC46425	Straight	#1
SC46426	Straight	#2
SC46427	Straight	#3
SC46428	Straight	#4
SC46429	Straight	#5
SC46430	Straight	#6

## COBB CURETTES, STAINLESS STEEL HANDLE

All stainless steel working end and stainless steel, knurled handle, 11" (28 cm).




CODE #	TIP SHAPE	SIZE
SC46433S	Angled	#3/0
SC46434S	Angled	#2/0
SC46435S	Angled	#0
SC46436S	Angled	#1
SC46437S	Angled	#2
SC46438S	Angled	#3
SC46439S	Angled	#4
SC46440S	Angled	#5
SC46441S	Angled	#6

CODE #	TIP SHAPE	SIZE
SC46422S	Straight	#3/0
SC46423S	Straight	#2/0
SC46424S	Straight	#0
SC46425S	Straight	#1
SC46426S	Straight	#2
SC46427S	Straight	#3
SC46428S	Straight	#4
SC46429S	Straight	#5
SC46430S	Straight	#6

# SPECIALTY BONE CURETTES


## EPSTEIN CURETTES

Reverse angle, hollow handle, 8" (20.5 cm).



CODE #	DESCRIPTION
SC49580	#000
SC49582	#00
SC49584	#0

## SCOVILLE CURETTES




CODE #	DESCRIPTION
SC49574	Straight, 10" (25.5 cm)
SC49576	Down angle, 10" (25.5 cm)
SC49578	Up angle, 10" (25.5 cm)

## JANSEN CURETTE



CODE #	DESCRIPTION
SC29632	Double-ended, oval cups, 6" (15 cm)

## MARTINI BONE CURETTE



CODE #	DESCRIPTION
SC29620	Double-ended, 4 mm and 5 mm round cups, 5-1/2" (14 cm)

## HATFIELD BONE CURETTES


Stainless steel with satin finish. Oval shaped tips, overall length 9" (23 cm). Complete set of nine curettes



CODE #	TIP SHAPE	SIZE	TIP LENGTH
584160	000	2.5 mm	4.0 mm
	Straight	2.8 mm	4.6 mm
	Straight	3.3 mm	5.0 mm
	Straight	3.8 mm	6.0 mm
	Straight	4.8 mm	7.6 mm
	Straight	5.8 mm	9.2 mm
	Straight	7.2 mm	11.6 mm
	Straight	8.7 mm	12.8 mm
	Straight	10.5 mm	15.5 mm

## VOLSCANN CURETTES

Double-ended, oval and round cups.



CODE #	DESCRIPTION
SC46560	5" (12.5 cm)
SC46561	8" (20.5 cm)
SC46562	6-1/2" (16.5 cm)


## WILLIGER BONE CURETTE



CODE #	DESCRIPTION
SC29630	Double-ended, 3 mm and 4 mm oval cups, 5-1/2" (14 cm)

## VOLSCANN CURETTES


Double-ended, oval cups.



CODE #	DESCRIPTION
SC49200	Double-ended, straight, 5" (12.5 cm)
SC49202	Double-ended, straight, 5 1/2" (14 cm)
SC46550	Double-ended, oval cups, 6 1/2" (16.5 cm)
SC49204	Double-ended, straight, 8" (20.5 cm)

## RING CURETTES


Hollow handle with fenestrated cup.



CODE #	DESCRIPTION
SC55680	Size 0, small, 10" (25.5 cm)
SC55684	Size 0, small, 11-1/2" (29 cm)
SC55681	Size 2, medium, 10" (25.5 cm)
SC55685	Size 2, medium, 11-1/2" (29 cm)
SC55682	Size 6, large, 10-1/4" (26 cm)
SC55686	Size 6, large, 12" (30.5 cm)

## CONE RING CURETTES


Knurled handle and fenestrated cup, 9" (23 cm).



CODE #	TIP SHAPE	SIZE	
SC49450	Straight	#1	8 mm diameter
SC49452	Straight	#2	6 mm diameter
SC49454	Straight	#3	3 mm diameter
SC49456	Angled 35 degrees	#1	8 mm diameter
SC49458	Angled 35 degrees	#2	6 mm diameter
SC49460	Angled 35 degrees	#3	3 mm diameter

# RAY PITUITARY CURETTES

7-1/2" (19 cm).





CODE #	TIP SHAPE	SIZE
SC49403	Vertical angle 45 degrees	3 mm
SC49404	Vertical angle 45 degrees	4 mm
SC49405	Vertical angle 45 degrees	5 mm
SC49413	Vertical angle 90 degrees	3 mm
SC49414	Vertical angle 90 degrees	4 mm
SC49415	Vertical angle 90 degrees	5 mm
SC49423	Horizontal angle 45 degrees	3 mm
SC49424	Horizontal angle 45 degrees	4 mm
SC49425	Horizontal angle 45 degrees	5 mm
SC49433	Horizontal angle 90 degrees	3 mm
SC49434	Horizontal angle 90 degrees	4 mm
SC49435	Horizontal angle 90 degrees	5 mm

# ENDOSCOPIC CURETTES


## ENDO BRUN CURETTES

Endoscopic with phenolic handle



CODE #	TIP WIDTH	TIP	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
SC55532	Angled 40°	2 x 4 mm cup	2.7"	9.3"	0.285"	5"	17" (43 cm)
SC55523	Straight	3 x 6 mm cup	2.7"	9.3"	0.285"	5"	17" (43 cm)
SC55533	Angled 40°	3 x 6 mm cup	2.7"	9.3"	0.285"	5"	17" (43 cm)
SC55543	Angled 40°	3 x 6 mm cup	3-1/2"	10"	0.45"	6-1/2"	20" (51 cm)
SC55556	Angled 30°	6 x 9 mm cup	3-1/2"	10"	0.45"	6-1/2"	20" (51 cm)


# ENDO EPSTEIN CURETTES



CODE #	TIP WIDTH	TIP	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
SC55593	Reverseangled 40°	3 x 6 mm	2.7"	9.3"	0.285"	5"	17" (43 cm)
SC55586	Reverseangled 30°	6 x 9 mm	3-1/2"	10"	0.45"	6-1/2"	20" (51 cm)

# ENDO CONE RING CURETTES


Endoscopic with phenolic handle



CODE #	TIP WIDTH	TIP	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
SC55573	Angled 40°	3 mm	2.7"	9.3"	0.285"	5"	17" (43 cm)
SC55566	Angled 30°	6 mm	10"	10"	0.45"	6-1/2"	20" (51 cm)

# ENDO SCOVILLE CURETTE

Endoscopic with phenolic handle.



CODE #	TIP WIDTH	TIP	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	OVERALL LENGTH
SC55510	Straight	9 x 15 mm	3-1/2"	10"	0.45"	6-1/2"	20" (51 cm)

# SMARTRACK CURETTE TRAYS

## BRUN SMARTRACK FOR 7" CURETTES

Aluminum and stainless steel.



CODE #	DESCRIPTION
SC98160	Holds 12 7" SCedic® Brun Cures

## COBB SMARTRACK FOR 11" INSTRUMENTS

Aluminum and stainless steel.



CODE #	DESCRIPTION
SC98170	Holds nine 11" SCedic® Cobb Cures, Elevators or Gouges

## HIBBS AND BRUN SMARTRACK FOR 9" INSTRUMENTS

Aluminum and stainless steel.



CODE #	DESCRIPTION
SC98180	Holds 12 9" SCedic® Brun Cures, or Hibbs Chisels, Gouges or Osteotomes