




Orthopedic Surgical Sets Catalogue 2026




6.0mm Spinal Pedicle Screw Set

Description:

- Awl
- Probe, Straight
- Probe, Curved
- Probe, Straight
- Double end Probe
- Fixed pin driver
- Ball fixed pin
- Cylinder fixed pin
- Tap
- Tap
- Tap
- Tap
- Ratchet Handle
- Multi-axial Screw Hex driver
- Fix angle screws driver
- T-handle
- Wrenching bar
- Plug Starter, Hex
- Pre-Plug Driver
- 20 cm Template rod
- Rod Bender
- Rod holder
- Sagittal Plane In-Situ Bender
- Sagittal Plane In-Situ Bender
- Rod Pusher
- Rod Reducer Plier
- Rocker
- Counter Reduction Sleeve Wrench
- Final Plug Driver
- SW 4.0 Driver shaft

- Curved Compressor
- Curved Distractor
- Rod Rotation Wrench
- Multi-axial Reduction Screw Hex driver
- Fix angle Reduction screws driver
- Rod Reducer
- Plier-Reduction Extension sleeve A
- Rod Reducer
- Plier-Reduction Extension sleeve B
- Reduction screw Rocker
- Reduction Tab Breaker
- Rod Rotation Forcep
- 10Nm Torsion Limited Driver
- Crosslink Meter
- Crosslink driver
- Crosslink Plate Bender
- Crosslink plate holder
- Crosslink Plate Locker
- Screw Driver
- Crosslink plate holder (Φ 3.0)
- Instrument case Box

SC-OS-101



6.5mm Cannulated Screw Instruments Set

Description:

- Depth Gauge (0–120 mm)
- Threaded Guider Wire 2.5 mm
- Threaded Guider Wire 2.5 mm
- Cannulated Drill Bit with Limited Block 4.5 mm
- Cannulated Countersink Φ 9
- Hex Key
- Wrench for Adjustable Parallel Wire Guider
- Multiple Wire Guider
- Tap Cannulated Screw 6.5 mm
- Screwdriver Hexagonal 3.5 mm
- Cleaning Stylet 2.5 mm
- Drill Sleeve
- Adjustable Parallel Wire Guider
- Cannulated Screwdriver Hexagonal 3.5 mm
- Aluminum Box




SC-OS-102

ACI-PCL Reconstruction Guide Set Of 5 Pcs Knee Orthopedic Instrument

Description:

- Material: A Grade Stainless Steel
- Rusting Prevention Procedure: Passivated
- Ultrasonic Cleaned: Yes
- Lubricate: Yes
- Polished: No (but possible)
- Sand-Blasted: No (but possible)
- Dull-Polished: Yes (German Dull)
- TC: No
- Gold Platted: HALF HANDLE
- Usage: Left Hand or Right Hand
- Tests Performed: Boil Test, Performance Test, Shape Test
- Packing: Individually Packed




SC-OS-103

ACL PCL Knee Arthroscopy Surgical Surgery Instruments Set 43 Pcs

Description:

- Main zig 1
- Elbow aimer 1
- PCL Tibial 1
- PCL Femoral 1
- Angle Buliate 1
- Graft passing 2.4mm Drill Bit 1
- Graft passing top (Beacht pin) 1
- Cannulated Tibial Reamer: Dia.6.0mm 1
- Cannulated Tibial Reamer: Dia.7.0mm 1
- Cannulated Tibial Reamer: Dia.8.0mm 1
- Cannulated Tibial Reamer: Dia.9.0mm 1
- Cannulated Tibial Reamer: Dia.10.0mm 1
- Cannulated Tibial Reamer: Dia.11.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.6.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.7.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.8.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.9.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.10.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.11.0mm 1
- Depth Gauge for ACL/PCL reconstruction 1
- Femoral Aimer, 5mm 1
- Femoral Aimer, 7mm 1
- Tunnel Plug, 7-9 mm 1
- Tunnel Plug, 8-9 mm 1
- Tunnel Plug, 9-10 mm 1
- Tunnel Dilator, 7mm 1
- Tunnel Dilator, 8mm 1
- T-Handle for Tunnel Dilators 1
- Cannulated screwdriver for Biol.cc interference screw 1
- Impactor for spiked ligament staple 1
- Arthroscopic Probe 1
- Microfracture Awl, 30° 1
- Microfracture Awl, 45° 1
- Microfracture Awl, 90° 1
- Graft Prep Station 1
- Graft Sizing Block 1
- Tendon Stripper, Closed end, 6mm 1
- Tendon Stripper, Open end, 7mm 1
- Straight Rasp 1
- Convex Rasp 1
- Femoral PCL Elevator 1

- Trocar Pin 1
- Empty Container for Knee Arthroscopy Instruments 1



SC-OS-104


ACL PCL Knee Arthroscopy Surgical Surgery Instruments Set 48 Pcs

Description:

- Main zig 1
- Elbow aimer 1
- PCL Tibial 1
- PCL Femoral 1
- Angle Bullater 1
- Graft passing 2.4mm Drill Bit 1
- Graft passing pin (Beastin pin) 1
- Cannulated Tibial Reamer: Dia.6.0mm 1
- Cannulated Tibial Reamer: Dia.7.0mm 1
- Cannulated Tibial Reamer: Dia.8.0mm 1
- Cannulated Tibial Reamer: Dia.9.0mm 1
- Cannulated Tibial Reamer: Dia.10.0mm 1
- Cannulated Tibial Reamer: Dia.11.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.6.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.7.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.8.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.9.0mm 1

- Cannulated Femoral Flowertip Reamer: Dia.10.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.11.0mm 1
- Depth Gauge for ACL/PCL reconstruction 1
- Femoral Aimer, 5mm 1
- Femoral Aimer, 7mm 1
- Tunnel Plug, 7-8 mm 1
- Tunnel Plug, 8-9 mm 1
- Tunnel Plug, 9-10 mm 1
- Tunnel Dilator, 7mm 1
- Tunnel Dilator, 8mm 1
- Tunnel Dilator, 9mm 1
- T-Handle for Tunnel Dilators 1
- Cannulated screwdriver for BioLoc interference screw 1
- Impactor for spiked ligament staple 1
- Arthroscopic Probe 1
- Microfracture Awl, 30° 1
- Microfracture Awl, 45° 1
- Microfracture Awl, 90° 1
- Graft Prep Station 1
- Graft Sizing Block 1
- Tendon Stripper, Closed end, 6mm 1
- Tendon Stripper, Open end, 7mm 1
- Straight Rasp 1
- Convex Rasp 1

- Femoral PCL Elevator 1
- Trocar Pin 1
- Empty Container for Knee Arthroscopy Instruments 1




SC-OS-105

Anterior Cervical Fixation System Instrument Set Orthopedic

Description:

- Plate and Screw Box
- Brace Screwdriver
- Drill Bit Ø2.5
- Screwdriver Hex SW2.5
- Drill Bit Ø2.5
- Hand Drill
- Screwdriver Sleeve
- Plate Holder Forcep
- Screwdriver Sleeve
- Cervical Hook
- Surette
- Guider
- Screwdriver Hex SW2.5
- Distractor
- Tap with Stopper
- Location Screw Screwdriver
- AWL
- Brace Screw
- Nut Screwdriver
- Plate Bender
- Location Screw
- Aluminium Box



SWI SC-OS-106

Anterior Cervical Plate System Instrument Set For Spine Surgery

Description:

- Plate and Screw Box
- Drill Bit Ø2.5
- Drill Bit Ø2.5
- Screwdriver Sleeve
- Screwdriver Sleeve
- Curette
- Screwdriver Hex SW2.5
- Tap with Stopper
- AWL
- Nut Screwdriver
- Brace Screwdriver
- Screwdriver Hex SW2.5
- Hand Drill
- Plate Holder Forcep
- Cervical Hook
- Cervical Hook
- Guider
- Distractor
- Location Screw Screwdriver
- Brace Screw
- Plate Bender
- Location Screw
- Aluminium Box



SWI SC-OS-107

SWIFT CUT INTERNATIONAL

Anterior Cervical Retractor Kit

Description:

- Retractor Blade Curved 20*30
- Retractor Blade Curved 20*30
- 4 Toothed Retractor Blade 20*30
- 4 Toothed Retractor Blade 20*30
- Retractor Blade Curved 24*40
- Retractor Blade Curved 24*40
- 5 Toothed Retractor Blade 24*40
- 5 Toothed Retractor Blade 24*40
- Retractor Blade Blunt 24*40
- Retractor Blade Blunt 24*40
- Retractor
- Retractor Blade Curved 20*50
- Retractor Blade Curved 20*50
- 4 Toothed Retractor Blade 24*40
- 4 Toothed Retractor Blade 24*40
- Retractor Blade Curved 25*60
- Retractor Blade Curved 25*60
- 5 Toothed Retractor Blade 25*60
- 5 Toothed Retractor Blade 25*60
- Retractor Blade Blunt 25*60
- Retractor Blade Blunt 25*60
- Aluminium Box



SWI SC-OS-108

Arthroscopy Shoulder Instruments Set

Description:

- SHOULDER INSTRUMENT CSCE
- PROBE, HOOK 5.4mm, TIP W/5mm MARKINGS
- BIRDBEAK, 2.75mm, 45° UP TIP
- BIRDBEAK, 2.75mm, 22° UP TIP
- SUTURE RETRIEVER, ø3.4mm, STRAIGHT
- SHOULDER DEBRIDEMENT RSCP
- KNOT PUSHER, CLOSED END
- SLAP RSCP
- GLENOID RSCP
- SHOULDER TISSUE ELEVATOR
- SHOULDER TISSUE ELEVATOR
- ROTATOR CUFF GRSCPER, 4.2mm
- CURETTE BOTH SIDE CUT, 5.4mm RING
- PENETRATOR SUTURE RETRIEVER, 15°
- EXTRA LONG SWITCHING STICK
- SUTURE HOOK
- SUTURE CUTTER, OPEN ENDED, LEFT NOTCH
- FIBER TAPE CUTTER
- FIBER TAPE RETRIEVER W/SR HANDLE
- PUNCH FOR 4.5mm PUSH LOCK & CORK SCREW
- CORK SCREW PUNCH
- 5.5/6.5MM SWIVELOCK
- PUNCH FOR 3mm PUSH LOCK
- SCORPION SUTURE PSCSER, HUMPBAC
- SPEAR, TROOFLITE OBITURATOR
- DRILL FOR SHORT 2.9mm PUSH LOCK
- DRILL FOR 3.5mm PUSH LOCK
- PassPort Hemostat, Curved
- PassPort Measuring Guide
- Aluminium Tray




SWIFTCUT SC-OS-109

SWIFTCUT
INTERNATIONAL

Austin Moore Set With Sterilize Case

Description:

- Austin moore guage 01
- Austin moore impactor 01
- Rasp narrow 01
- Rasp broad 01
- Tommy bars 02
- Bone skid 01
- Box Chisel 01
- Cork screw 01
- Nail extractor set (contain 2 weights, 1 rod , 2 hooks) 05 stainless steel. Total Quantity. With BOX 14 Pcs




SWIFTCUT INTERNATIONAL SC-OS-110

Black Belt ACD Cervical Retractor Set With Titanium Blades

Description:

- 1 Retractor Frame RIGHT
- 1 Retractor Frame LEFT
- 2 Blade Insertion Clamps
- 14 Promged Blades 34mm-70mm
- 14 Blunt Blades 34mm-70mm
- 12 Two Level Promged Blades 40mm-55mm
- 6 Two Level Blunt Blades 40mm-55mm
- 1 Container



SC-OS-111

Bone Curette Set Orthopedic Surgical Instrument With Box

Description:

- Purple 4/0--- Straight
- Purple 4/0--- Angled
- Red 2/0--- Straight
- Red 2/0--- Angled
- Blue 0--- Straight
- Blue 0--- Angled
- Gold 2--- Straight
- Gold 2--- Angled



SC-OS-112

Bookwalter Kit Retractor System Surgical Premium Grade Quality Set With Box

Description:

- Ring + Blade TRAY
- 1"POST COUPLING BALFOUR RETRACTOR
- HORIZONTAL FLEXBAR
- HARRINGTON FLEXBAR
- MEDIUM OVAL RING
- SEGMENTEDRING
- MALLEABLE RETRACTOR 2 (2"x6"/51MMX152MM)
- MALLEABLE RETRACTOR 1 (3"x6"/76MMX152MM)
- KELLYRETRACTOR 2 (2"x3"/51MMX76MM)
- KELLYRETRACTOR 1 (2"x4"/51MMX110MM)
- KELLYRETRACTOR 1 (2"x5"/51MMX127MM)
- KELLYRETRACTOR 1 (2"x6"/51MMX152MM)
- 6 TILT RATCHET MECHANISMS
- 2 RATCHET MACHANISMS




SC-OS-113

Bookwalter Kit Retractor System Surgical Premium Grade Quality Set With Box

Description:

- Ring + Blade TRAY
- 1" POST COUPLING BALFOUR RETRACTOR
- HORIZONTAL FLEXBAR
- HARRINGTON FLEX BAR
- MEDIUM OVAL RING
- SEGMENTED RING
- GELPI RETRACTOR
- MALLEABLE RETRACTOR 2 (2"x6"/51MMX152MM)
- MALLEABLE RETRACTOR 1 (3"x6"/76MMX152MM)
- KELLYRETRACTOR 2 (2"x3"/51MMX76MM)
- KELLYRETRACTOR 1 (2"x4"/51MMX110MM)
- KELLYRETRACTOR 1 (2"x5"/51MMX127MM)
- KELLYRETRACTOR 1 (2"x6"/51MMX152MM)
- 6 TILT RATCHET MECHANISMS
- 2 RATCHET MACHANISMS





SC-OS-114

Broken Screw Removal Instrument Set

Description:

- Status Screwdriver, 6mm
- Stable Screwdriver, 6mm
- Spump Screwdriver, 1.5mm
- Tipping Screwdriver, 1.5mm
- Trepin, 4mm - Trepin, 6mm
- Toggle, 6.5mm
- Hexagon Screwdriver, 1.5mm
- Hexagon Screwdriver, 2.0mm
- Hexagon Screwdriver, 3.0mm
- Hexagon Screwdriver, 5.5mm
- Hexagon Screwdriver, 4.0mm
- Hexagon Screwdriver, 5.5mm
- Hexagon Screwdriver, 6.0mm
- Hexagon Screwdriver, 5.5mm
- Trepin, 8.2mm - Trepin, 8.2mm
- Trepin, 6.0mm - Trepin, 8.0mm
- Trepin, 6.5mm - Trepin, 6mm
- Hexagon Slip Extractor, 1.2mm
- Hexagon Slip Extractor, 2.5mm
- Hexagon Slip Extractor, 2.7mm
- Hexagon Slip Extractor, 3.5mm

- Hexagon Slip Extractor, 4.0mm
- Hexagon Slip Extractor, 5.5mm
- Torx Screwdriver, T18
- Torx Screwdriver, T19
- Torx Screwdriver, T20
- Torx Screwdriver, T21
- Torx Screwdriver, T22
- Torx Screwdriver, T23
- Torx Screwdriver, T24
- Sohn / Style Signum
- Drill Bit, 1mm - Drill Bit, 5mm - Drill Bit, 8
- Resurer - Sharp Hook + Quick Handle, 1
- Quick Handle, Straight
- Trepin, 4.5mm - Trepin, 3.5mm
- Bone Screw Holding Fotops
- Bone Screw
- Sterilization Case Overall Dimensions: 15 inch L x 5.8 inch W x 9.28 inch H




SC-OS-115

Caspar Lumbar Micro Discectomy Retractor Set

Description:

- 1 Sterilization Container
- 1 Casper Directory
- 1 Blade Ejector Forceps
- 7 Lateral Blades - 40 mm, 45 mm, 50 mm, 55 mm, 60 mm, 65 mm, 70 mm
- 7 Speculas - 40 mm, 45 mm, 50 mm, 55 mm, 60 mm, 65 mm, 70 mm




SC-OS-116

DHS&DCS Orthopedic Instrument Set Aluminum Wit Box

Description:

- DHS/DCS Reamer 1PC
- DHS/DCS Reamer 1PC
- DHS/DCS Tap 12.5mm 1PC
- DHS/DCS Angle/guide 150' 1PC
- DHS/DCS Angle/guide 45' 1PC
- DHS/DCS Impector 1PC
- Depth guage Large 1PC
- Screw Driver 4.5mm 1PC
- Drill Sleeve 3.5mm 4.5mm 1PC
- Drill Sleeve 4.5mm 6.5mm 1PC
- Drill Guide 3.5mm 4.5mm 1PC
- Missling guage 1PC
- DHS/DCS Pana 1PC
- DHS/DCS Pana sleeve 2PC
- Drill bit 3.5mm 2PC
- Drill bit 4.5mm 2PC
- Tap Handle 4.5mm
- Tap Handle 4.5mm



SC-OS-117

Distal Radius 2.4 mm Orthopedic Best Instrument Set Aluminum Box Screw And Plate

Description:

- Sterilization box
- 10 Plates with (2.7 mm forking holes) Material 316 L (Stainless Steel)
- Locking screws for for 10 Plates. Material 316 L (Stainless Steel)
- Screw Driver
- 2 Bone liver
- 2 Sleeves
- 3 Drill bits
- Depth gauge
- Sleeves handle
- Drill & tap guider


Total Number of Screw (110 Screws)

Cortical Screws: 316 L (Stainless Steel)

- 8mm (5), 10mm (5), 12mm (5), 14mm (5), 16mm (5), 18mm (5), 20mm (5), 22mm (5), 24mm (5), 26mm (5), 28mm (5)

Locking Screws: 316 L (Stainless Steel)

- 8mm (5), 10mm (5), 12mm (5), 14mm (5), 16mm (5), 18mm (5), 20mm (5), 22mm (5), 24mm (5), 26mm (5), 28mm (5)




SC-OS-118

Femoral Interlocking Intramedullary Nails Set For Femur Bone Fracture Surgery

Description:

- Universal Joint 1
- Drill Sieve 2
- Open Wrench1
- Quick Coupling T-Handle 1
- Screwdriver 1
- Drill Sieve 1
- Hex Wrench1
- Hex Wrench1
- Drill Sieve 1
- Drill Sieve1
- T-handle Drill Bit 2
- Depth Gauge 10-100mm 1
- Drill Bit 5.2*300mm 1
- Drill Bit $\phi 4.0 \times 300$ 1
- Drill Bit with Limitator $\phi 4.0 \times 300$ 1
- Reamer 10mm 1
- Reamer 11mm 1
- Reamer 12mm 1
- Reamer 13mm 1
- AWL 1
- Guide Wire 1
- Hammer 1
- Guide Rod 1
- Connector Rod 1
- Nail Connector 1
- Location Rod 1
- Sighting Device Handle1

- Bolt 1
- Bolt 1
- Bolt 1
- Bolt 1
- Fixed Connector 1
- Location Forcep 1
- Connector1
- Distal Guider 1
- Box 1



SC-OS-119

Flexarm Nathen Liver Retractor Set





SC-OS-120

Humerus Interlocking Nail Set With Box

Description:

- Ring + Blade TRAY
- 1" POST COUPLING BALFOUR RETRACTOR
- HORIZONTAL FLEX BAR
- HARRINGTON FLEX BAR
- MEDIUM OVAL RING
- SEGMENTED RING
- GELPI RETRACTOR
- MALLEABLE RETRACTOR 2(2"x6"/51MMX152MM)
- MALLEABLE RETRACTOR 1(3"x6"/76MMX152MM)
- KELLY RETRACTOR 2(2"x3"/51MMX76MM)
- KELLY RETRACTOR 1(2"x4"/51MMX110MM)
- KELLY RETRACTOR 1(2"x5"/51MMX127MM)
- KELLY RETRACTOR 1(2"x6"/51MMX152MM)
- 6 TILT RATCHET MECHANISMS
- 2 RATCHET MACHANISMS



SWIFT SC-OS-121

SWIFT CUT
INTERNATIONAL

SWIFT CUT
INTERNATIONAL

SWIFT CUT
INTERNATIONAL

SWIFT CUT
INTERNATIONAL



SWIFT CUT
INTERNATIONAL

SWIFT CUT
INTERNATIONAL

SWIFT CUT
INTERNATIONAL

SWIFT CUT
INTERNATIONAL

SWIFT CUT
INTERNATIONAL




SWIFT CUT
INTERNATIONAL

K-Wire Instruments Complete Set With Sterilize Case

Description:

- 1Wire Tightener Handle with two pegs1
- 2Plier cum cutter1
- 3Wire bending pliers1
- 4Wire Passer2
- 5Vice Grip1
- 6Forceps for holding cerclage wire2
- 7Wire with eye of 1.0mm dia., length 280mm10
- 8Wire with eye of 1.25mm dia., length 280mm10
- 9Wire coil, 1 each of 1.0 and 1.25mm dia.2
- 10K Wire Dia. 1.0mm, Length 150mm10
- 11K Wire Dia. 1.25mm, Length 150mm10
- 12K Wire Dia. 1.5mm, Length 150mm10
- 13K Wire Dia. 2.0mm, Length 150mm10
- 14K Wire Dia. 2.5mm, Length 150mm10
- 15K Wire Dia. 3.0mm, Length 150mm10




SC-OS-122

Large Fregment Instruments Complete Set With Sterilize Case

Description:

- Drii | Bit $\Phi 3.2 \times 150\text{mm}$ 2
- Drii | Bit $\Phi 4.5 \times 150\text{mm}$ 2
- Countersink $\Phi 8$ 1
- Tap canvas 4.5mm 1
- Tap cancelions 6.5mm 1
- Screw Forceps 1
- 4.5 mm DCP Drill Guide neutral and load 1
- Double Drill Sleeve $\Phi 3.2/6.5$ 1
- Universal Drill Guide $\Phi 3.2/4.5$ 1
- Template 11x12u0 1
- Template 12x155 1
- Template 12x212 1
- Depth Gauged-90 1
- Extraction Screw For Hexagonal 3.5mm, conical 1
- Hollow Reamer 8mm 1
- Screwdiver Hexagonal 3.5mm 1
- Reduction Forceps with points, broad, ratchet, 210mm 1
- Reduction Forceps with serrated Jaw, ratchet, 210mm 1
- Self-centering Bone Holding Forceps.260mm 2
- Periosteal Elevator 8mm, round 1
- Periosteal Elevator 16mm, flat 1
- Holmium Retractor, Width 44mm 2
- Bending Iron L/R 1
- quick coupling t handle 1
- screwdriver shaft 4.5l
- Surgical Sterilization Box




SC-OS-123

Leyla Brain Retractor Set

Description:

- 2 Socket Joints
- 1 Coupling Head
- 1 Holding Rod
- 2 Flexible Arms with fix Base
- 2 Flat Spatulas with Holders



SC-OS-124

MARTIN ARM Neurosurgery Retractor Liver System Complete Orthopedic

Description:

- Martin Arm Surgical Retractor System Complete Set
- BRAND NEW Quality A+
- 1 Bed Support
- 1 Rod
- 1 Moveable Arm
- 3 Wires as shown in Picture



SC-OS-125

Micro Mini Instruments Complete Set With Sterilize Case

Description:

DRILL BITS & GUIDES:

- Drill bit: 1.1mm (Qty: 2)
- Drill bit: 1.5mm (Qty: 2)
- Drill bit: 2.0mm (Qty: 2)
- Drill bit: 2.5mm (Qty: 2)
- Drill bit: 3.2mm (Qty: 2)

Neutral & Loaded Drill Guide (Qty: 1 each):

- 1.1mm & 1.5mm
- 1.5mm & 2.0mm
- 2.0mm & 2.7mm
- 2.5mm & 3.5mm
- 3.2mm & 4.0mm

TAP SLEEVES & TAPS:

- Tap sleeve: 1.1mm & 1.5mm (Qty: 1)
- Tap sleeve: 2.7mm & 2.0mm (Qty: 1)
- Tap sleeve: 2.5mm & 3.5mm (Qty: 1)
- Tap sleeve: 3.2mm & 4.0mm (Qty: 1)
- HA 1.5 Tap (Qty: 1)
- HA 2.0 Tap (Qty: 1)
- HA 2.7 Tap (Qty: 1)
- HA 3.5 Tap (Qty: 1)
- HA 4.0 Tap (Qty: 1)

MESCURING TOOLS:

- Depth Gauge (Qty: 1)
- Counter sink large (Qty: 1)
- Counter sink small (Qty: 1)

SCREWDRIVERS & TOOLS:


- Hex Head Screw driver 1.5mm (Qty: 1)
- Hex head Screw driver 2.7mm (Qty: 1)
- Hex head screw driver 3.5mm (Qty: 1)
- Plate Benders pair-modern design (Qty: 1)
- Plate benders small (Qty: 1)
- Periosteal Elevator (Round Head) (Qty: 1)
- Broken screw externaliser (Qty: 1)
- Periosteal Elevator (Flat Head) (Qty: 1)
- Reduction Forceps with point (Qty: 1)
- Reduction Forceps (Qty: 1)
- Self-centering forceps (Qty: 1)

STORAGE & IMPLANTS:

- Case (Qty: 1)
- Screw rack-3.5/4.0mm (Qty: 1)
- Screw rack-1.5/2.0/2.7 (Qty: 1)
- Bone plate 1.5mm (Qty: 1)
- Bone plate 2.0mm (Qty: 1)
- Bone plate 2.7mm (Qty: 1)
- Bone plate 3.5/4.0mm (Qty: 1)

Cortical screws self-tapping (Qty: 1 each):

- 1.5mm
- 2.0mm
- 2.7mm
- 3.5mm
- 4.0mm-cancellous



SC-OS-126

Mis Tublar Retractor Spine system Complete Set

Description:

- 16 mm Cups (Qty: 7)
- 18 mm Cups (Qty: 7)
- 20 mm Cups (Qty: 7)
- Dilators (5 mm × 50 mm) (Qty: 7)
- Complete Bed Support System with Malleable Arm (Qty: 1)




SC-OS-127

Mis Tublar Retractor Spine System Complete Set With Sterilize Case

Description:

- (7) 16 mm Cups
- (7) 18 mm Cups
- (7) 20 mm Cups
- (7) Dilators (5 mm x 50 mm)(1)
- Complete Bed Support System with Malleable Arm




SWIFT SC-OS-128

SWIFTCUT INTERNATIONAL

Mueller-Type Cement Removal Instruments Set

Description:

- Narrow Cement Removal Gouge, Short – Shaft Length: 150mm – Gouge: 9mm, negative (Qty: 1)
- Narrow Cement Removal Gouge, Long – Shaft Length: 240mm – Gouge: 9mm, negative (Qty: 1)
- Narrow Offset Cement Removal Gouge – Shaft Length: 240mm – Gouge: 9mm, negative (Qty: 1)
- Acetabular Chisel – Shaft Length: 240mm – Chisel: 7.5mm (Qty: 1)
- Offset Chisel – Shaft Length: 150mm – Chisel: 9mm (Qty: 1)
- Flared Angle Gouge – Shaft Length: 240mm – Gouge: 9mm, positive, angle 15° downward (Qty: 1)
- Wide Gouge – Shaft Length: 240mm – Gouge: 11.5mm, negative (Qty: 1)
- "V" Splitter – V-Shaped Chisel: 7mm (Qty: 1)
- Cement Splitting Osteotome – Shaft Length: 240mm (Qty: 1)
- Saddle Punch – Shaft Length: 240mm – Punch: 16.5mm x 6.5mm (Qty: 1)
- Cement Removal Osteotome, Short – Shaft Length: 150mm – Osteotome: 8mm (Qty: 1)
- Cement Removal Osteotome, Long – Shaft Length: 240mm – Osteotome: 8mm (Qty: 1)
- 4.4mm Drill – 4.4mm (Qty: 1)



SC-OS-129

Omni Track Complete Retractor Set

Description:

- Main zig 1
- Elbow aimer 1
- PCL Tibial 1
- PCL Femoral 1
- Angle Bullater 1
- Graft passing 2.4mm Drill Bit 1
- Graft passing pin (Beastin pin) 1
- Cannulated Tibial Reamer: Dia.6.0mm 1
- Cannulated Tibial Reamer: Dia.7.0mm 1
- Cannulated Tibial Reamer: Dia.8.0mm 1
- Cannulated Tibial Reamer: Dia.9.0mm 1
- Cannulated Tibial Reamer: Dia.10.0mm 1
- Cannulated Tibial Reamer: Dia.11.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.6.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.7.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.8.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.9.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.10.0mm 1
- Cannulated Femoral Flowertip Reamer: Dia.11.0mm 1
- Depth Gauge for ACL/PCL reconstruction 1
- Femoral Aimer, 5mm 1
- Femoral Aimer, 7mm 1
- Tunnel Plug, 7-8 mm 1
- Tunnel Plug, 8-9 mm 1
- Tunnel Plug, 9-10 mm 1
- Tunnel Dilator, 7mm 1
- Tunnel Dilator, 8mm 1
- Tunnel Dilator, 9mm 1
- T-Handle for Tunnel Dilators 1
- Cannulated screwdriver for BioLoc interference screw 1
- Impactor for spiked ligament staple 1
- Arthroscopic Probe 1
- Microfracture Awl, 30° 1
- Microfracture Awl, 45° 1
- Microfracture Awl, 90° 1
- Graft Prep Station 1
- Graft Sizing Block 1
- Tendon Stripper, Closed end, 6mm 1
- Tendon Stripper, Open end, 7mm 1
- Straight Rasp 1
- Convex Rasp 1



SC-OS-130

Orthopedic Bone Screw Driver Quick Coupling Handle 9 Pieces Set

Description:

- Screw Driver 1.5 mm (Qty: 1)
- Screw Driver 2.0 mm (Qty: 1)
- Screw Driver 2.5 mm (Qty: 1)
- Screw Driver 3.0 mm (Qty: 1)
- Screw Driver 3.5 mm (Qty: 1)
- Screw Driver 4.0 mm (Qty: 1)
- Screw Driver 4.5 mm (Qty: 1)
- Screw Driver 5.0 mm (Qty: 1)
- Includes: Quick Coupling handle (Qty: 1)



SC-OS-131

Orthopedic Hudson Brace Hand Drill Set of 7 Pcs

Description:

- 1 Brace. Length 10.2°
- 1 Extension piece for brace. Length 3.9°
- 1 Cushing 14mm Flat Drill. Length 4°
- 1 Burr, diameter 9mm, Cylindrical. Length 4°
- 1 Burr, diameter 14mm, Cylindrical. Length 3.9°
- 1 Burr, diameter 16mm, Spherical. Length 3.9°
- 1 Burr, diameter 22mm, Spherical. Length 3.9°
- Premium Quality Box



SC-OS-132

Screws Cortical 4.5mm Locking 5.0mm Orthopedics Screw With Box

Description:

- Cortical screw 4.5 mm Size each 18 mm To 80 mm 156
- Locking screw 5.0 mm Size each 18 mm TO 90 mm 150
- Cancellous Screw 6.5 mm Half Thread Size each 30mm To 95 mm 44
- Cancellous screw 6.5mm Full Thread Size each 30mm To 95mm 44
- Orthopedic screws Lowt Of 394 peices Box Large Fragments
- Cortical 4.5mm Locking 5.0mm Cancellous 6.5mm 394PCS



SC-OS-133

Small Fregment Instruments Complete Set With Sterilize Case

Description:

- Elevator 9mm (Qty: 1)
- Elevator 12mm (Qty: 1)
- Bending Spanner (Qty: 2)
- Screw Extractor (Qty: 1)
- Double Drill Guide 2.5-3.5 (Qty: 1)
- Double Drill Guide 3.0-4.0 (Qty: 1)
- Double Drill Guide 2.5 (Qty: 1)
- Drill bit 2.5mm (Qty: 2)
- Drill bit 3.0mm (Qty: 2)
- Screw driver Hex SW2.5 (Qty: 1)
- Screw Tap HA3.5 (Qty: 1)
- Screw Tap HB4.0 (Qty: 1)
- Depth Gauges 0-60mm (Qty: 1)
- Countersink Drill (Qty: 1)
- Reduction Forceps (small) (Qty: 1)
- 1° Point Reduction Forceps (small) (Qty: 1)
- Bone Holding Forceps Self-Centering (Qty: 2)
- Antslip Screw Driver (Qty: 1)
- Screw Rack (Qty: 1)
- Aluminum Case (Qty: 1)



SC-OS-134

Spine Sugery Instruments Set With Aluminum Box

Description:

- Awl
- -Probe, Straight
- -Probe, Curved
- Probe, Straight
- Double end Probe
- fixed pin driver
- ball fixed pin
- cylinder fixed pin
- Tap
- Tap
- Tap
- Tap
- Ratchet Handle
- Multi-axial Screw Hex driver
- Fix angle screws driver
- T-handle
- Wrenching bar
- Plug Starter, Hex
- Pre-Plug Driver
- 20cm Template rod
- Rod Bender
- Rod holder
- Sagittal Plane In-Situ Bender
- Sagittal Plane In-Situ Bender
- Rod Pusher
- Rod Reducer Piler
- Rocker
- Counter Reduction Sleeve Wrench
- Final Plug Driver
- SW4.0 Driver shaft
- Curved Compressor

- Curved Distractor
- Rod Rotation Wrench
- Multi-axial Reduction Screw Hex driver
- Fix angle Reduction screws driver
- Rod Reducer Piler- Reduction Extension sleeve A
- Rod Reducer Piler- Reduction Extension sleeve B
- Reduction screw Rocker
- Reduction Tab Breaker
- Rod Rotation Forcep
- 10Nm Torsion Limited Driver
- Crosslink Meter
- Crosslink driver
- Crosslink Plate Bender
- Crosslink plate holder
- Crosslink Plate Locker Screw Driver
- Crosslink plate holder($\Phi 3.0$)
- Instrument case Box



SC-OS-135

Thompson Retractor Complete Set Orthopedic Surgical Instrument Stainless Steel

Description:

- Rail Clamps 26"(55cm) tall with 2 cam joint - 2 PCS
- Crossbars Angled arm 18"(8" * 10") - 3 PCS
- Bilateral crossbars 7 1/4" * 13 1/2" * 7 1/4" - 1 PCS
- Bilateral Crossbar Hinged 11" * 8 1/2" * 11" - 1 PCS
- Retractor Handles S-Lock (non-angling) 8" - 3 PCS
- Retractor Handles Lock Clip-on Angling 10" - 1 PCS
- General Surgery Blades 76mm * 51mm - 1 PCS
- General Surgery Blades 83mm * 76mm - 2 PCS
- Kelly Blades 64mm * 76mm - 1 PCS
- Kelly Blades 76mm * 89mm - 1 PCS
- Malleable Blades 13mm * 90mm - 1 PCS
- Malleable Blades 13mm * 152mm - 1 PCS
- Malleable Blades 13mm * 203mm - 1 PCS
- Malleable Blades 25mm * 203mm - 1 PCS
- Malleable Blades 25mm * 104mm - 1 PCS
- Malleable Blades 51mm * 203mm - 1 PCS
- Malleable Blades 64mm * 254mm - 1 PCS
- Malleable Blades 76mm * 254mm - 1 PCS
- Malleable Blades 102mm * 254mm - 1 PCS
- St.Marks Blades 64mm * 178mm - 1 PCS
- Harrington Blades 64mm * 152mm - 1 PCS
- Wide Balfour Blades 83mm*57mm - 1 PCS
- Wide Balfour Blades 102mm*57mm - 1 PCS
- Wide Balfour Blades 114mm*57mm - 1 PCS



SC-OS-136



DEDICATED
SALES
TEAM

CONTACT US NOW!

FOR EXCEPTIONAL CUSTOMER CARE.

DETAILS ✓
QUOTATION ✓
VIRTUAL MEETING ✓
ORDER ✓

CONTACT US NOW!

Location: 
Shahabpura, Sialkot, Pakistan

Website: 
www.swiftcutint.com

Email: 
info@swiftcutint.com

Phone Number: 
+92 305 91976121