

OSTEOTOMES



Osteotomes are used to cut or reshape bone with controlled force during orthopedic surgeries. Designed with sharp, tapered blades, they allow for precise bone sectioning while minimizing trauma to surrounding tissue. Made from high-strength surgical steel, our osteotomes offer excellent durability, balance, and grip for safe and effective use in a variety of procedures.



LAMBOTTE


MINI-LAMBOTTE OSTEOTOMES, STRAIGHT 5"

Straight, without calibration lines, 5" (12.5 cm).



CODE #	DESCRIPTION
SC46336	1/8" (4 mm)
SC46338	1/4" (6 mm)
SC46340	5/16" (8 mm)
SC46342	3/8" (10 mm)
SC46344	1/2" (12 mm)
SC46343	9/16" (14 mm)
SC46345	5/8" (15 mm)
SC46346	11/16" (18 mm)
SC46347	3/4" (19 mm)
SC46335	13/16" (20 mm)

LAMBOTTE OSTEOTOMES, STRAIGHT 7"



CODE #	DESCRIPTION
SC46348	1/8" (4 mm)
SC46349	3/16" (5 mm)
SC46350	1/4" (6 mm)
SC46358	5/16" (8 mm)
SC46351	3/8" (10 mm)
SC46352	1/2" (12 mm)
SC46353	5/8" (16 mm)
SC46354	3/4" (19 mm)
SC46355	7/8" (22 mm)
SC46356	1" (25 mm)

LAMBOTTE-STYLE OSTEOTOMES


flat thin blade, 5" (12.5 cm)



CODE #	DESCRIPTION
051082	2 mm
051084	4 mm
051086	6 mm
051088	8 mm
051090	10 mm
051092	12 mm
051095	15 mm


LAMBOTTE OSTEOTOMES, STRAIGHT 9"

Straight, without calibration lines, 9" (23 cm).



CODE #	DESCRIPTION
SC46304	1/4" (6 mm)
SC46306	1/2" (13 mm)
SC46310	3/4" (19 mm)
SC46314	1" (25 mm)
SC46318	1-1/4" (32 mm)
SC46322	1-1/2" (38 mm)


LONG BEVEL OSTEOTOMES



CODE #	DESCRIPTION
SC69940	Straight, 10 mm
SC69946	Straight, 20 mm
SC69950	Straight, 25 mm
SC69970	Curved, 5 mm
SC69975	Curved, 10 mm

LAMBOTTE OSTEOTOMES, CURVED 9"

Curved, without calibration lines, 9" (23 cm).



CODE #	DESCRIPTION
SC46328	1/4" (6 mm)
SC46329	1/2" (12 mm)
SC46330	3/4" (19 mm)
SC46331	1" (25 mm)
SC46332	1-1/4" (32 mm)
SC46333	1-1/2" (38 mm)

SPECIALTY OSTEOTOMES

Flat thin blade, 6" (15 cm).



CODE #	DESCRIPTION
051102	2 mm
051103	3 mm
051104	4 mm
051105	5 mm
051106	6 mm
051107	7 mm
051110	10 mm

MINI OSTEOTOMES



Curved, 5" (12.5 cm).



CODE #	DESCRIPTION
441207	4 mm
441208	6 mm

HOKE HOKE OSTEOTOMES

5-1/4" (13.5 cm).




CODE #	TIP	TIP WIDTH
SC46497	Straight	1/8" (3 mm)
SC46498	Straight	3/16" (4 mm)
SC46499	Straight	1/4" (6 mm)
SC46500	Straight	5/16" (8 mm)
SC46501	Straight	3/8" (10 mm)
SC46502	Straight	1/2" (13 mm)
SC46503	Straight	5/8" (16 mm)
SC46504	Straight	3/4" (19 mm)

CODE #	TIP	TIP WIDTH
SC46700	Curved	1/8" (3 mm)
SC46701	Curved	3/16" (4 mm)
SC46702	Curved	1/4" (6 mm)
SC46703	Curved	5/16" (8 mm)
SC46704	Curved	3/8" (10 mm)
SC46706	Curved	1/2" (13 mm)
SC46708	Curved	5/8" (16 mm)
SC46709	Curved	3/4" (19 mm)

SILVER

SILVER OSTEOTOMES

Single guard, 7" (18 cm)



CODE #	DESCRIPTION
SC58100	Straight
SC58102	Guard, right
SC58104	Guard, left

U.S. Army

U.S. ARMY PATTERN OSTEOTOMES

7" (18 cm)




CODE #	DESCRIPTION
SC46284	6 mm
SC46286	12 mm
SC46288	18 mm

SHEEHAN

SHEEHAN OSTEOTOMES

Straight, 6-1/4" (16 cm)




CODE #	DESCRIPTION
SC46268	2 mm
SC46270	4 mm
SC46272	6 mm
SC46274	8 mm
SC46276	10 mm
SC46278	12 mm
SC46280	15 mm

STILLE

STILLE OSTEOTOMES

Straight, 8" (20.5 cm).



CODE #	DESCRIPTION
SC49002	10 mm
SC49004	14 mm
SC49006	18 mm
SC49008	22 mm

HIBBS


HIBBS OSTEOTOMES



CODE #	DESCRIPTION
SC1011	

HIBBS ENDOSCOPIC OSTEOTOME

Straight tip, 18-1/2" (47 cm).




CODE #	TIP WIDTH	TIP LENGTH	SHAFT LENGTH	SHAFT DIAMETER	HANDLE LENGTH	HANDLE LENGTH
SC55642	3/8" (10 mm)	3-1/2"	10"	0.43" (11 mm)	10-3/4"	4-1/4"

SMITH PETERSON

SMITH PETERSON OSTEOTOME


8" (20.5 cm).



CODE #	tip	TIP WIDTH
SC46448	Straight	1/4"
SC46449	Straight	3/8"
SC46450	Straight	1/2"
SC46451	Straight	5/8"
SC46452	Straight	3/4"
SC46453	Straight	1"
SC46454	Straight	1-1/4"
SC46457	Curved	1/4"
SC46458	Curved	3/8"
SC46459	Curved	1/2"
SC46460	Curved	5/8"
SC46461	Curved	3/4"
SC46462	Curved	1"
SC46463	Curved	1-1/4"

INTERCHANGEABLE OSTEOTOME, CHISEL AND GOUGE SET

INTERCHANGEABLE OSTEOTOME, CHISEL AND GOUGE SET




CODE #	DESCRIPTION
SC49630	Complete set with handle, key and all blades
SC49632	Phenolic handpiece only, 7-3/8" (18.5 cm)
SC49634	Key only
SC49668	Osteotome blade, 5 mm
SC49672	Osteotome blade, 10 mm
SC49676	Osteotome blade, 16 mm
SC49680	Osteotome blade, 25 mm
SC49652	Chisel blade, 5 mm
SC49656	Chisel blade, 10 mm
SC49640	Gouge blade, 1" wide, 50 mm radius
SC49644	Gouge blade, 1" wide, 60 mm radius
SC49648	Gouge blade, 1" wide, 70 mm radius

COBB

COBB OSTEOTOME

All stainless steel, 11" (28 cm).



CODE #	tip	TIP WIDTH
SC49350	Straight	1/4"
SC49351	Straight	3/8"
SC49352	Straight	1/2"
SC49353	Straight	5/8"
SC49354	Straight	3/4"
SC49355	Straight	1"
SC49356	Straight	1-1/4"
SC49357	Curved	1/4"
SC49358	Curved	3/8"
SC49359	Curved	1/2"
SC49360	Curved	5/8"
SC49361	Curved	3/4"
SC49362	Curved	1"
SC49363	Curved	1-1/4"