

# RONGEURS

Rongeurs are used for trimming, shaping, and removing small pieces of bone or soft tissue in orthopedic procedures. Featuring sharp, scoop-shaped jaws, they provide precise cutting action with minimal effort. Made from high-quality surgical stainless steel, our Rongeurs offer smooth operation, durability, and a firm grip for controlled and efficient performance in delicate surgical areas.



# SINGLE ACTION DESIGN


## SYNOVECTOMY RONGEURS

2 x 5 mm bite, extremely delicate




CODE #	DESCRIPTION
SC47252	Slightly curved, 6-7/8" (17.5 cm)
SC47253	Strongly curved, 6-3/4" (17 cm)

## RONGEUR #4A



CODE #	DESCRIPTION
SC47275	Curved, 4 x 7 mm bite, 6-1/2" (16.5 cm)

## RONGEUR #4



CODE #	DESCRIPTION
SC47274	Curved, 4 x 9 mm bite, 5-1/2" (14 cm)

## MICRO FRIEDMAN RONGEURS


1.5 x 5 mm bite. For soft tissue and cancellous bone only. Extremely delicate.



CODE #	DESCRIPTION
SC47240	Straight, 5-3/4" (14.5 cm)
SC47242	30 degree angled, 5-1/2" (14 cm)
SC47244	90 degree angled, 5-1/2" (14 cm)

## MINI FREIDMAN FONGEURS

Delicate, 2 x 8 mm bite, 6" (15 cm)



CODE #	DESCRIPTION
SC47263	Straight
SC47264	Curved

## KLEINERT-KUTZ SYNOVECTOMY RONGEUR



CODE #	DESCRIPTION
065202	25 degrees strongly curved, 5-1/4" (13.5 cm)

## FRIEDMAN RONGEUR



CODE #	DESCRIPTION
SC47262	Curved, delicate, 4 x 11 mm bite, 5-3/4" (14.5 cm)

## BANE RONGEUR

Curved, 3 x 16 mm bite, 7" (18 cm)




CODE #	DESCRIPTION
SC47310	25 degrees strongly curved, 5-1/4" (13.5 cm)

## DEAN RONGEUR




CODE #	DESCRIPTION
SC47286	Curved, 3 mm bite, 5-3/4" (14.5 cm)

## BACON RONGEUR



CODE #	DESCRIPTION
SC47308	Curved on side, 5 x 15 mm bite, 8" (20.5 cm)

## CHICHIRELLI (BLUMENTHAL) RONGEUR



CODE #	DESCRIPTION
SC47267	Curved 45 degrees, 4 x 13 mm bite, 6" (15 cm)

## LEMPERT RONGEURS

Delicate, 3 x 7 mm bite




CODE #	DESCRIPTION
SC47258	Straight, 6" (15 cm)
SC47260	Curved, 6" (15 cm)
SC47254	Straight, 7-1/2" (19 cm)
SC47256	Curved, 7-1/2" (19 cm)

## BEYER BONE BITING RONGEUR



CODE #	DESCRIPTION
028120	Curved jaw, 16 x 3 mm cup, 6-3/4" (17 cm)

## MEAD RONGEUR



CODE #	DESCRIPTION
SC47276	Curved, 5 x 13 mm bite, 6-1/2" (16.5 cm)

## HARTMANN RONGEUR



CODE #	DESCRIPTION
SC47285	Curved, 4 x 14 mm bite, 7-1/4" (18.5 cm)

## ADSON RONGEURS

7 x 15 mm bite, 8" (20.5 cm).



CODE #	DESCRIPTION
SC47302	Straight
SC47306	Curved

## DEAN (BLUMENTHAL) RONGEURS

4 x 11 mm, 6" (15 cm).



CODE #	DESCRIPTION
SC47266	Curved 30 degrees
SC47268	Curved 90 degrees

## LEUR RONGEURS

8 x 12 mm bite, 7-1/8" (18 cm).



CODE #	DESCRIPTION
SC47278	Straight
SC47284	Curved

## ADSON RAT TOOTH RONGEUR

Teeth facilitate easier bone excision.



CODE #	DESCRIPTION
SC47285	Straight with teeth, 8 x 15 mm bite, 8" (20.5 cm)

## DOUBLE ACTION DESIGN

## KLEINERT-KUTZ SYNOVECTOMY RONGEUR



CODE #	DESCRIPTION
SC47553	Curved on side, delicate, 2 x 10 mm bite, 6" (15 cm)

## KLEINERT-KUTZ SYNOVECTOMY RONGEURS

Delicate, 6-1/8" (15.5 cm).




CODE #	DESCRIPTION
SC47560	Curved, 2 mm bite
SC47562	Curved, 3 mm bite

## KLEINERT-KUTZ SYNOVECTOMY RONGEUR



CODE #	DESCRIPTION
SC47558	Fully curved, delicate, 3 mm bite, 6-1/8" (15.5 cm)

## KLEINERT-KUTZ DOUBLE ACTION RONGEUR



CODE #	DESCRIPTION
065212	10 degrees lightly curved, 18.5 x 5 mm cup, 7-1/4" (18.5 cm)

## BEYER RONGEURS

Curved, 7" (18 cm).



CODE #	DESCRIPTION
SC47520	3 x 14 mm bite
SC47522	4 x 15 mm bite

## RUSKIN RONGEURS

6" (15.5 cm).



CODE #	DESCRIPTION
SC47523	Straight, 2 x 14 mm bite
SC47532	Straight, 3 x 15 mm bite
SC47533	Curved, 2 x 15 mm bite
SC47534	Curved, 3 x 15 mm bite

## RUSKIN RONGEUR

7-1/4" (18.5 cm).



CODE #	DESCRIPTION
SC47526	Straight, 4 x 16 mm bite
SC47528	Straight, 6 x 16 mm bite
SC47531	Curved, 6 x 16 mm bite

## RUSKIN RONGEUR



CODE #	DESCRIPTION
SC47525	Curved, 5 x 10 mm bite, 7-1/4" (18.5 cm)

## RUSKIN RONGEUR



CODE #	DESCRIPTION
341836	Double action, curved, 4 mm bite, 6-1/8" (15.5 cm)

## ZAUFAL-JANSEN RONGEUR



CODE #	DESCRIPTION
SC47524	Curved, 5 x 15 mm bite, 7" (18 cm)

## SMITH-PETERSON LAMINECTOMY RONGEURS

3 x 16 mm bite, 9-1/2" (24 cm).



CODE #	DESCRIPTION
SC47474	Straight
SC47476	Curved
SC47478	Full curved

## SMITH-PETERSON LAMINECTOMY RONGEUR



CODE #	DESCRIPTION
SC47525	Double action, curved 15 degrees, 7-1/2" (19 cm)

## STILLE-LUER ECHLIN DUCKBILL RONGEURS

Curved on side, 9" (23 cm)



CODE #	DESCRIPTION
SC47472	6 x 13 mm bite
SC47473	6 x 15 mm bite

## LEKSELL RONGEURS


Curved 50 degrees, 9-1/2" (24 cm).



CODE #	DESCRIPTION
SC47494	3 x 13 mm bite
SC47495	4 x 13 mm bite
SC47497	5 x 13 mm bite
SC47498	6 x 13 mm bite
SC47499	7 x 13 mm bite
SC47496	8 x 16 mm bite


## LEKSELL RONGEURS

Curved 15 degrees, 9" (23 cm)



CODE #	DESCRIPTION
SC47492	3 x 16 mm bite
SC47486	4.5 x 16 mm bite
SC47487	5 x 16 mm bite
SC47488	6 x 16 mm bite
SC47489	7 x 16 mm bite
SC47490	8 x 16 mm bite

## STILLE-LEUR RONGEURS



CODE #	DESCRIPTION
SC47458	Straight, 10 x 18 mm, 8-1/2" (21.5 cm)
SC47462	Straight, 10 x 18 mm, 10-1/2" (26.5 cm)
SC47460	Curved, 10 x 18 mm, 8-3/8" (21 cm)
SC47464	Curved, 10 x 18 mm, 10-1/2" (26.5 cm)
SC47461	Curved on side, 10 x 16 mm, 8-1/2" (21.5 cm)

# ECHLIN DUCKBILL RONGEURS


Curved on side, 9" (23 cm)



CODE #	DESCRIPTION
SC47480	2 x 12 mm bite
SC47482	3 x 12 mm bite
SC47484	4 x 12 mm bite

# SYPERT RONGEURS


14" (35.5 cm).



CODE #	DESCRIPTION
SC47380	4 x 16 mm
SC47382	5 x 16 mm
SC47384	8 x 16 mm
SC47386	12 x 16 mm


# SPINAL RONGEURS AND PUNCHES

## HANDLE STYLES



SHAFT LENGTH measured in closed position


## IDENTIFICATION CHART



# CLEANPUNCH RONGEURS

- The first Kerrison that can be completely cleaned—eliminating all debris between shafts
- Precise cutting action— surgery proceeds smoothly
- Easily disassembles and reassembles— all parts remain connected

**Patent Number: 5,961,531**



## HANDLING THE CLEANPUNCH RONGEUR

### TO OPEN:

1. Hold punch in left hand. Turn knob to “open” position (align arrows).
2. Push knob up with left thumb. Grip top slider on serrations with right hand and pull back.
3. Clean interior surfaces.

### TO CLOSE:

1. Put knob in “open” position.
2. Align guides on top slider.
3. Push top slider forward until it clicks into place.
4. Turn knob to “close” position.
5. Test for proper functioning

## FEATURES AND BENEFITS:

Patented design enables the instrument to be completely cleaned, eliminating the “buildup” of bone, tissue and fluids that cause sticking and resistance to action.

Provides a consistently smooth cutting action.

Eliminates distractions that interrupt a surgeon’s concentration required in delicate spinal procedures.

Instant identification enables surgeons to quickly locate the right bite size. Large, legible mm sizes are etched on handles, as well as color coded.

## COLOR-CODED FOR INSTANT BITE SIZE RECOGNITION

- 1 mm – Red
- 2 mm – Blue
- 3 mm – Yellow
- 4 mm – Black
- 5 mm – Green
- 6 mm – White

# KERRISON RONGEURS — CLEANPUNCH 2

## KERRISON CLEANPUNCH 2 LAMINECTOMY RONGEURS

Handle style VII, thin footplate, 7" (18 cm)



CODE #	BITE	SIZE
SC59300	Forward	1 mm
SC59303	Forward	2 mm
SC59306	Forward	3 mm
SC59309	Forward	4 mm
SC59312	Forward	5 mm

## KERRISON CLEANPUNCH 2 LAMINECTOMY RONGEURS

Handle style VII, thin footplate, 8" (20.5 cm).



CODE #	BITE	SIZE
SC59313	Forward	1 mm
SC59315	Forward	2 mm
SC59316	Forward	3 mm
SC59321	Forward	4 mm
SC59322	Forward	5 mm

## CLEANPUNCH CLEANING SYSTEM

Holds six CleanPunches in open position for maximum exposure during cleaning.


Lightweight anodized aluminum construction.

Minimize footprint – compact design has removable lid that fits under rack.

When open, instrument handles are exposed for quick identification and easy access in the OR.

Silicone-coated brackets prevent metal to metal contact. Protects sharp cutting edge during cleaning.

Fits inside a standard container or can be wrapped separately



CODE #	DESCRIPTION
SC47525	Approximately 13-1/2" x 7" x 6"

# KERRISON RONGEURS – THIN FOOTPLATE DESIGN

## CERVICAL KERRISON THIN FOOTPLATE RONGEURS

Thin footplate designed specifically for small bone in the cervical area.


Gold hexagonal screw denotes thin footplate.

Bite size etched on handle for quick identification.

Care should be taken when using delicate thin footplate punches.

Do not twist or pry instrument while cutting

Do not use on large bones.



See the Instructions For Use for all instructions and contraindications.

Handle style VII, thin footplate.


CODE #	BITE	SIZE	OVERALL LENGTH
SC55891	Forward	1 mm	6" (15 cm)
SC55892	Forward	2 mm	6" (15 cm)
SC55895	Forward	3 mm	6" (15 cm)
SC55898	Forward	4 mm	6" (15 cm)
SC55901	Forward	5 mm	6" (15 cm)
SC55904	Forward	1mm	7" (18 cm)
SC55908	Forward	2mm	7" (18 cm)
SC55911	Forward	3mm	7" (18 cm)
SC55914	Forward	4mm	7" (18 cm)
SC55917	Forward	5mm	7" (18 cm)
SC55925	Forward	1mm	8" (20.5 cm)
SC55926	Forward	2 mm	8" (20.5 cm)
SC55927	Forward	3 mm	8" (20.5 cm)
SC55928	Forward	4 mm	8" (20.5 cm)
SC55929	Forward	5mm	8" (20.5 cm)

CODE #	BITE	SIZE	OVERALL LENGTH
KM55905	Up	1 mm	7" (18 cm)
KM55909	Up	2 mm	7" (18 cm)
KM55912	Up	3 mm	7" (18 cm)
KM55915	Up	4 mm	7" (18 cm)
KM55918	Up	5 mm	7" (18 cm)
KM55906	Down	1mm	7" (18 cm)
KM55910	Down	2mm	7" (18 cm)
KM55913	Down	3mm	7" (18 cm)
KM55916	Down	4mm	7" (18 cm)
KM55919	Down	5mm	7" (18 cm)

# KERRISON RONGEURS – STANDARD DESIGN

## LOVE-KERRISON LAMINECTOMY PUNCHES

Handle style IV, large handle, regular footplate



CODE #	BITE	SIZE	OVERALL LENGTH
SC47676	Forward	2 mm	8" (20.5 cm)
SC47670	Forward	3 mm	6" (15 cm)
SC47671	Forward	3 mm	7" (18 cm)
SC47678	Forward	3 mm	8" (20.5 cm)
SC47680	Forward	4 mm	8" (20.5 cm)
SC47674	Forward	5 mm	6" (15 cm)
SC47682	Forward	5 mm	8" (20.5 cm)
SC47689	Forward	6 mm	11" (28 cm)

## LOVE-KERRISON RONGEURS


Handle style II, large handle, regular footplate, 8" (20.5 cm).



CODE #	BITE	SIZE
SC47620	Forward	1 mm
SC47621	Forward	2 mm
SC47622	Forward	3 mm
SC47623	Forward	4 mm
SC47624	Forward	5 mm

## LOVE-KERRISON LAMINECTOMY PUNCHES


Handle style IV, large handle, regular footplate.



CODE #	BITE	SIZE	OVERALL LENGTH
SC47644	Up	3 mm	6" (15 cm)
SC47645	Up	4 mm	6" (15 cm)
SC47646	Up	5 mm	6" (15 cm)
SC47648	Down	3 mm	6" (15 cm)
SC47649	Down	4 mm	6" (15 cm)
SC47650	Down	5 mm	6" (15 cm)

# LOVE-KERRISON LAMINECTOMY PUNCHES


Handle style IV, large handle, regular footplate.



CODE #	BITE	SIZE	OVERALL LENGTH
KM47644	Up	3 mm	6" (15 cm)
KM47645	Up	4 mm	6" (15 cm)
KM47646	Up	5 mm	6" (15 cm)
KM47648	Down	3 mm	6" (15 cm)
KM47649	Down	4 mm	6" (15 cm)
KM47650	Down	5 mm	6" (15 cm)
KM47652	Up	3 mm	8" (20.5 cm)
KM47853	Up	4 mm	8" (20.5 cm)
KM47654	Up	5 mm	8" (20.5 cm)
KM47656	Down	3 mm	8" (20.5 cm)
KM47657	Down	4 mm	8" (20.5 cm)
KM47658	Down	5 mm	8" (20.5 cm)

# SPURLING-KERRISON PUNCHES

Handle style IV, regular footplate, 8" (20.5 cm).



CODE #	BITE	SIZE
SC55870	Forward	3 mm
SC55871	Forward	4 mm
SC55872	Forward	5 mm

# SCHLESINGER CERVICAL LAMINECTOMY RONGEURS


Handle style IV, large handle, regular footplate, 6" (15 cm).



CODE #	BITE	SIZE
SC47610	Up	3 mm
SC47612	Up	5 mm
SC47614	Down	3 mm
SC47616	Down	5 mm


# FERRIS-SMITH KERRISON RONGEURS

Handle style IV, thin footplate, 7" (18 cm).



# FERRIS-SMITH KERRISON RONGEURS

Handle style V, regular footplate



# FERRIS-SMITH KERRISON PUNCHES

Ring handle (handle style VI), regular footplate.



# FERRIS-SMITH ENDOSCOPIC KERRISON RONGEURS

N RONGEURS Graduated tip, ring handle, 3-1/4" tip length, .470" shaft diameter, working length 14" (35.5 cm), overall 20" (51 cm).



# KERRISON RONGEUR



CODE #	BITE	SIZE	OVERALL LENGTH
411625	Up	3 mm	3-1/2" (9 cm)

# MOORELAND RONGEUR



CODE #	DESCRIPTION
SC47850	5 x 15 mm, 13-3/4" (35 cm)

# KERRISON MULTI-BITE PUNCHES

Fenestrated upper shaft ejects cut pieces of bone and tissue. Delicate, thin footplate, for cervical use only, 7" (18 cm).



CODE #	BITE	SIZE
SC55980	Up, fenestrated	2 mm
SC55981	Up, fenestrated	3 mm

# FORAMINOTOMY PUNCHES

Handle style V, regular footplate, 8" (20.5 cm).



CODE #	BITE	SIZE
SC48030	Curved up	2 mm
SC48033	Curved up	3 mm
SC48032	Curved up	4 mm

# IVD RONGEURS

## DECKER IVD RONGEUR



CODE #	BITE	SIZE	OVERALL LENGTH
SC55858	Straight	2 x 6 mm	6" (15 cm)

## SELVERSTONE IVD RONGEUR



CODE #	BITE	SIZE	OVERALL LENGTH
SC55860	Straight	2 x 6 mm	6" (15 cm)

## PEAPOD IVD RONGEUR



CODE #	BITE	SIZE	OVERALL LENGTH
SC55842	Up angle	2 x 6 mm	5-1/2" (14 cm)

## WILDE IVD RONGEURS

Fenestrated jaws, straight



CODE #	DESCRIPTION	OVERALL LENGTH
SC47707	Small, 4 x 9 mm bite	5" (12.5 cm)
SC47708	Medium, 5 x 10 mm bite	5" (12.5 cm)
SC47709	Large, 6 x 11 mm bite	5" (12.5 cm)

## CUSHING IVD RONGEURS



CODE #	BITE	SIZE	OVERALL LENGTH
SC47710	Straight	2 x 10 mm	5" (12.5 cm)
SC47730	Straight	2 x 10 mm	7" (18 cm)
SC47712	Forward	2 x 10 mm	5" (12.5 cm)
SC47732	Forward	2 x 10 mm	7" (18 cm)
SC47714	Down angle	2 x 10 mm	5" (12.5 cm)
SC47734	Down angle	2 x 10 mm	7" (18 cm)

# SPURLING IVD RONGEURS



CODE #	BITE	SIZE	OVERALL LENGTH
SC47800	Straight	4 x 10 mm	5" (12.5 cm)
SC47840	Straight	4 x 10 mm	7" (18 cm)
SC47802	Forward	4 x 10 mm	5" (12.5 cm)
SC47842	Forward	4 x 10 mm	7" (18 cm)
SC47804	Down angle	4 x 10 mm	5" (12.5 cm)
SC47844	Down angle	4 x 10 mm	7" (18 cm)

# CASPAR IVD RONGEURS



CODE #	BITE	SIZE	OVERALL LENGTH
SC55935	Forward	2 x 12 mm	6" (15 cm)
SC55936	Forward	3 x 12 mm	6" (15 cm)
SC55937	Forward	4 x 14 mm	6" (15 cm)
SC55938	Forward	5 x 14 mm	6" (15 cm)
SC55960	Forward	2 x 12 mm	6" (15 cm)
SC55961	Forward	3 x 12 mm	7" (18 cm)
SC55962	Forward	4 x 14 mm	7" (18 cm)
SC55963	Forward	5 x 14 mm	7" (18 cm)
SC55964	Forward	6 x 16 mm	7" (18 cm)
SC55972	Forward	3 x 12 mm	7" (18 cm)
SC55973	Forward	4 x 14 mm	8" (20.5 cm)
SC55974	Forward	5 x 14 mm	8" (20.5 cm)
SC55966	Forward	3 x 12 mm	8" (20.5 cm)
SC55967	Forward	4 x 14 mm	8" (20.5 cm)

# LOVE-GRUENWALD IVD RONGEURS


Fenestrated jaw, serrated cup



CODE #	BITE	SIZE	OVERALL LENGTH
SC47901	Straight	3 x 10 mm	5" (12.5 cm)
SC47902	Straight	3 x 10 mm	7" (18 cm)
SC47903	Forward	3 x 10 mm	5" (12.5 cm)
SC47904	Forward	3 x 10 mm	7" (18 cm)
SC47905	Down angle	3 x 10 mm	5" (12.5 cm)
SC47906	Down angle	3 x 10 mm	7" (18 cm)

# CASPAR IVD RONGEURS

6" (15 cm).



CODE #	BITE	SIZE
SC55946	Straight	2 x 10 mm
SC55947	Straight	3 x 10 mm
SC55948	Straight	4 x 10 mm
SC55949	Straight	5 x 10 mm
SC55952	Up angle	2 x 10 mm
SC55953	Up angle	3 x 10 mm
SC55954	Up angle	4 x 10 mm
SC55940	Down angle	2 x 10 mm
SC55941	Down angle	3 x 10 mm
SC55942	Down angle	4 x 10 mm

# SCHLESINGER IVD RONGEURS

9" (23 cm).



CODE #	BITE	SIZE
SC55690	Straight	3 x 10 mm
SC55691	Straight	5 x 10 mm
SC55693	Up angle	3 x 10 mm
SC55694	Up angle	5 x 10 mm

# SCHLESINGER IVD RONGEURS

Serrated jaws.



CODE #	BITE	SIZE	OVERALL LENGTH
SC47700	Straight	2 x 10 mm	6" (15 cm)
SC47702	Straight	3 x 10 mm	6" (15 cm)
SC47704	Straight	4 x 10 mm	6" (15 cm)
SC47703	Up angle	3 x 10 mm	6" (15 cm)

# OLDBERG IVD RONGEUR



CODE #	BITE	SIZE	OVERALL LENGTH
SC55856	Round	6 mm	7" (18 cm)

# FERRIS-SMITH IVD RONGEURS


Ring handle lessens chance of slippage. Handle style VI, 7" (18 cm) shaft.



CODE #	BITE	SIZE
SC47858	Straight	2 x 10 mm
SC47864	Straight	3 x 10 mm
SC47870	Straight	4 x 10 mm
SC47878	Straight	6 x 12 mm
SC47860	Up angle	2 x 10 mm
SC47866	Up angle	3 x 10 mm
SC47872	Up angle	4 x 10 mm
SC47873	Up angle	6 x 12 mm
SC47862	Down angle	2 x 10 mm
SC47868	Down angle	3 x 10 mm
SC47874	Down angle	4 x 10 mm
SC47875	Down angle	6 x 12 mm

# FERRIS-SMITH ENDOSCOPIC IVD RONGEURS

Handle style VI, ring handle, 3-1/4" tip length, 0.470" (12 mm) shaft diameter, 10-3/4" shaft length, 6" handle length, 20" (51 cm).



CODE #	BITE	SIZE
SC55762	Straight	2.5 x 10 mm
SC55763	Straight	4 x 10 mm
SC55767	Up angle	4 x 10 mm
SC55777	Straight, serrated	8 x 12 mm