




Our orthopedic screwdrivers are designed for precise control and secure handling of all standard implant screws. With ergonomic handles and torque-optimized shafts, they allow for accurate insertion and removal without stripping or slippage. Available in a variety of tip designs and sizes, each screwdriver ensures compatibility with a wide range of orthopedic systems.


SCREWDRIVERS

HEX HEAD HEXAGONAL SCREWDRIVERS



CODE #	DESCRIPTION
SC47952	2.5 mm screwdriver with sleeve, 7-7/8" (20 cm)
SC48368	2.5 mm screwdriver with notch, 8-1/4" (21 cm)
SC48349	3.5 mm screwdriver with notch, 9-5/8" (24.5 cm)
SC47955	Screw holding sleeve only for SC47952, 2" (5 cm)
SC48365	Screw holding sleeve only for SC48368, 3" (8 cm)
SC47983	Screw holding sleeve only for SC48349, 5" (12.5 cm)
SC47956	Screwdriver only for SC47952

LARGE HEXAGONAL SCREWDRIVER (PELVIC)



CODE #	DESCRIPTION
SC48353	3.5 mm screwdriver with notch, 11-3/4" (30 cm)
SC48364	Screwdriver only, 2.5 mm with notch, 11" (27 cm)
SC47983	Sleeve only for SC48353, 5" (12.5 cm)

HEXAGONAL SCREWDRIVER SHAFT (PELVIC)



CODE #	DESCRIPTION
SC47951	2.5 mm with notch, quick coupling end, 6-1/2" (16.5 cm)
SC47953	2.5 mm without notch, quick coupling end, 4" (10 cm)
SC47961	3.5 mm with notch, quick coupling end, 6-1/2" (16.5 cm)
SC48365	Sleeve only for SC47951, 3" (8 cm)
SC47983	Sleeve only for SC47961, 5" (12.5 cm)

SCREWDRIVER SLEEVES



CODE #	DESCRIPTION
SC48365	Sleeve only for SC47951, SC48348, SC48364, SC48368, 3" (8 cm)
SC47983	Sleeve only for SC48349, SC48353 and SC47961, 5" (12.5 cm)

LARGE HEXAGONAL SCREWDRIVER



CODE #	DESCRIPTION
SC48252 3	3.5 mm, quick coupling end, 3-7/8" (10 cm)

LARGE HEX SCREWDRIVER WITH UNIVERSAL JOINT



CODE #	DESCRIPTION
SC48396	3.5 mm tip, 11-5/8" (29.5 cm)

HEXAGONAL SCREWDRIVER SHAFT




CODE #	DESCRIPTION
SC48348	2.5 mm with notch, quick coupling end, 5-1/2" (14 cm)
SC48365	Sleeve only

LARGE HEXAGONAL SCREWDRIVER



CODE #	DESCRIPTION
SC47954	3.5 mm, quick coupling end, 3-7/8" (10 cm)

SCREWDRIVER




CODE #	DESCRIPTION
SC48670	1.5 mm hexagonal, 7-7/16" (200 mm)
SC48350	2.5 mm hexagonal, 9-27/32" (250 mm)
SC48352	3.5 mm hexagonal, 9-27/32" (250 mm)
SC48354	Single slot, 9-27/32" (250 mm)
SC48356	Cruciform, 9-27/32" (250 mm)
SC48358	Phillips, 9-27/32" (250 mm)

MINI HEXAGONAL SCREWDRIVER




CODE #	DESCRIPTION
SC47950	1.5 mm, mini quick coupling end, sleeve included, 2-1/4" (5.5 cm)
SC47950B	Shaft only, 1.5 mm with notch, mini quick coupling end, 2" (5 cm)
SC47959	Screw holding sleeve only

SINGLE SLOT SCREWDRIVER



CODE #	DESCRIPTION
SC48355	Phenolic handle, 6-3/4" (17 cm)

CRUCIFORM HEAD CRUCIFORM SCREWDRIVER



CODE #	DESCRIPTION
SC47957	1.5 / 2.0 mm screwdriver with sleeve, 7-1/2" (19 cm)
SC47959	Screw holding sleeve only for SC47957
SC47958	Screwdriver only for SC47957

MINI CRUCIFORM SCREWDRIVER BIT



CODE #	DESCRIPTION
SC45004	1.5 / 2 mm, mini quick coupling end, 2-1/2" (6.5 cm)

T-HANDLE LARGE HEX SCREWDRIVER W/ T-HANDLE

CODE #	DESCRIPTION
SC48395	3.5 mm, quick coupling end, 6-7/8" (17.5 cm)

DEPTH GAUGE SCREW DEPTH GAUGES

Used to determine screw length required during internal fracture fixation (mm graduated). Snaps apart for easy cleaning. IMPORTANT NOTE: When not in use, retract measuring end inside gauge for probe protection



CODE #	DESCRIPTION
SC47934	Mini, up to 30 mm, 4-3/4" (12 cm)
SC47936	Small, up to 50 mm, 6-1/4" (16 cm)
SC47940	Large, up to 100 mm, 8-3/4" (22 cm)
SC47935	Extra-long, up to 120 mm, 11" (28 cm)

INTERCHANGEABLE SCREWDRIVER SET

Screwdriver Sets



CODE #	DESCRIPTION
SC48473	Screwdriver handle. Allows K-Wire and pins to be inserted with Jacobs Universal Chuck SC48478. Cannulation 1/4" (6.3 mm), includes key, 5-1/2" (14 cm)
SC48474	Single slot, 5-1/2" (14 cm)
SC48475	Cruciform, 5-1/2" (14 cm)
SC48476	Phillips, 5-1/2" (14 cm)
SC48477	Small Phillips, 5-1/2" (14 cm)
SC48478	Jacobs Universal Chuck. Allows K-Wire and pins to be inserted with Jacobs Chuck. Cannulation 1/4" (6.3 mm), includes key, 3" (7.5 cm)
SC48478P	Key for SC48478 Jacobs Universal Chuck
SC48479	3.5 mm diameter, 5-1/2" (14 cm)
SC48480	2.5 mm diameter, 5-1/2" (14 cm)

SCREW HOLDING FORCEPS

SELF-RETAINING SCREW HOLDING FORCEPS



CODE #	DESCRIPTION
SC47930	3-1/2" (9 cm)

SELF-RETAINING SCREW HOLDING FORCEPS



CODE #	DESCRIPTION
SC45208	For 1.5 mm to 2.7 mm screws
SC45210	For 3.5 mm to 6.5 mm screws