




A comprehensive range of instruments for diagnostic and surgical procedures of the urinary tract and male reproductive system. Includes cystoscopes, resectoscopes, urethral dilators, stone-grasping forceps, biopsy forceps, catheters, and retractors. Designed for precise access, clear visualization, and atraumatic tissue handling within delicate, confined anatomy. Fine craftsmanship, polished stainless steel surfaces, and ergonomic grips ensure patient safety, surgeon control, and reliable long-term performance in both endoscopic and open urological surgery.

UROLOGY

Clutton Urethral Sounds

French gauge tapering, 265mm long




SC-UR/007	Set of 12 contains SC-UR/008 to SC-UR/008k
-----------	--

SC-UR/008	6 - 10	SC-UR/008f	18 - 22
SC-UR/008a	8 - 12	SC-UR/008g	20 - 24
SC-UR/008b	10 - 14	SC-UR/008h	22 - 26
SC-UR/008c	12 - 16	SC-UR/008i	24 - 28
SC-UR/008d	14 - 18	SC-UR/008j	26 - 30
SC-UR/008e	16 - 20	SC-UR/008k	28 - 32

Lister Urethral Sounds


English gauge, straight, 280mm long



SC-UR/010	0 - 1	SC-UR/010f	4 - 7
SC-UR/010a	0 - 2	SC-UR/010g	5 - 8
SC-UR/010b	0.5 - 3	SC-UR/010h	6 - 9
SC-UR/010c	1 - 4	SC-UR/010i	7 - 10
SC-UR/010d	2 - 5	SC-UR/010j	8 - 11
SC-UR/010e	3 - 6	SC-UR/010k	9 - 12

Winsbury White Urethral Bougie


French gauge, child size, 175mm long



SC-UR/063	1-4 to 3-6
SC-UR/064	5-8 to 7-10
SC-UR/065	9-12 to 11-14

Canny Ryall Meatal Urethral Bougie

English gauge, double-ended, tapering 160mm long



SC-UR/075	1-4 and 3-6
SC-UR/076	5-8 and 7-10
SC-UR/077	9-12 and 11-14
SC-UR/078	13-16 and 15-18
SC-UR/079	17-20 and 19-22
SC-UR/080	21-23 and 22-24
SC-UR/081	26-28 and 29-31

Buckston Browne Urethral Bougie

English gauge, 295mm long


SC-UR/020	4-6	SC-UR/027	11-13
SC-UR/021	5-7	SC-UR/028	12-14
SC-UR/022	6-8	SC-UR/029	13-15
SC-UR/023	7-9	SC-UR/030	14-16
SC-UR/024	8-10	SC-UR/031	15-17
SC-UR/025	9-11	SC-UR/032	16-18
SC-UR/026	10-12		



Clutton Urethral Bougie

French gauge, 260mm long


SC-UR/033	6 - 10	SC-UR/039	18 - 22
SC-UR/034	8 - 12	SC-UR/040	20 - 24
SC-UR/035	10 - 14	SC-UR/041	22 - 26
SC-UR/036	12 - 16	SC-UR/042	24 - 28
SC-UR/037	14 - 18	SC-UR/043	26 - 30
SC-UR/038	16 - 20	SC-UR/044	28 - 32



Guyon Urethral Bougie

French gauge, 275mm long


SC-UR/046	6-24
-----------	------



Lister Urethral Bougie

English gauge, 270mm long


SC-UR/047	4-7	SC-UR/053	10-13
SC-UR/048	5-8	SC-UR/054	11-14
SC-UR/049	6-9	SC-UR/055	12-15
SC-UR/050	7-10	SC-UR/056	13-16
SC-UR/051	8-11	SC-UR/057	14-17
SC-UR/052	9-12	SC-UR/058	15-18



Wyndham Powell Urethral Bougie


French gauge, 170mm long

SC-UR/059	12	SC-UR/062a	19	SC-UR/065	26
SC-UR/060	13	SC-UR/059b	20	SC-UR/066a	27
SC-UR/061	14	SC-UR/060b	21	SC-UR/067a	28
SC-UR/062	15	SC-UR/061b	22	SC-UR/068a	29
SC-UR/059a	16	SC-UR/062b	23	SC-UR/069a	30
SC-UR/060a	17	SC-UR/063	24	SC-UR/070a	31
SC-UR/061a	18	SC-UR/064	25	SC-UR/071a	32




Walton Bladder Retractor

SC-UR/001	Parallel action, with convex lateral blades - 25mm x 70mm and curved central blade - 70mm x 215mm
-----------	---



Thompson Walker Bladder Retractor

SC-UR/002	One posterior blade 95mm x 30mm, one anterior wall blade 175mm x 25mm tapering to 20mm, two pairs of lateral swivel blades, shallow: 55mm x 22mm, deep: 75mm x 25mm
-----------	---




Millin Retropubic Retractor

Selfretaining

SC-UR/003	
-----------	--


Millin retractor blades

SC-UR/004a	Short 25 x 55mm
SC-UR/004b	Long 25 x 80mm
SC-UR/004c	Centre 115 x 445mm



Millin Vulsellum Forceps

Screw joint for capsule, 225mm long



SC-UR/018	8-9 teeth, jaws 10mm wide
-----------	---------------------------

Millin Dissecting Spring Forceps

230mm long



SC-UR/014	3-4 teeth
-----------	-----------